

DETERMINATION OF ASH CONTENT

The ash amount is determined by burning off the plant and calculating the amount of ash obtained.

Amount of ash is determined;

- ❑ To understand if the drug is well purified,
- ❑ To determine the amount of inorganic substances, mineral salts, metals
- ❑ To make sure whether it is adulterated with other drugs.

- The value found is compared with the value specified for the drug in the pharmacopoeia, thus it is determined whether the drug is appropriate for the pharmacopoeia.
- Remaining part of a sample as ash is composed of inorganic substances. The organic part contains 'C' and is released as CO_2 during the burning process.

Constant weight of crucible: A crucible is weighed after it is left at a certain temperature (e.g. at a temperature of 600-800 °C in furnace) for a certain period of time until the difference between the last 2 weights is found to be ± 0.3 mg.

Experimental Procedure:

- Take out the crucibles, which are brought to constant weight in the furnace at 600°C , using a heated tongs.
- Allow to cool in a desiccator and weighed on precision balance. (P)
- Weigh 1.00 g of powdered sample in the crucible. (P1)
- Heat the crucible in a slightly inclined position on the triangular porcelain (for the homogeneous spread of heat and not to lose the ash) on the burner, burn off the sample over a light fire until the content of crucible is black and over a high heat until it turns into white ash. (pre-ashing process).
- Put the crucible in a muffle furnace at 600°C by a heated tongs, leave for 1 hour in the furnace, then allow the crucible at constant weight to cooldown in desiccator and weigh. (P2)

- Amount of sample : $P1 - P = A \text{ g}$
- Amount of ash: $P2 - P = B \text{ g}$
- $A \text{ g sample}$ $B \text{ g ash}$
- 100 $X = \% \text{ ash (a/a)}$