

AIM: to get familiar with the Linearly Modulated OSL (LM-OSL) signal and check the behavior of the LM-OSL measurements versus various stimulation parameters.

**Materials:** Aluminum oxide (TLD 500,  $\text{Al}_2\text{O}_3:\text{C}$ ) or Berilium oxide (BeO) or quartz

#### PROTOCOL:

Step 1. Irradiation (1 Gy)

Step 2. TL measurement (350 °C for synthetic materials, 500 °C for natural materials, HR=1 °C/s)

Step 3. Irradiation (1 Gy)

Step 4. LM- OSL measurement (Room temperature, 0% - 90% intensity, 200 s)

Step 5. TL measurement (350 °C for synthetic materials, 500 °C for natural materials, HR=1 °C/s) after OSL.

Step 6. Repeat steps 3-5 for various measurement temperatures (40 °C, 50 °C, 60 °C, 70 °C, 80 °C, 90 °C and 100 °C).

Step 7. Repeat steps 3-5 for various stimulation times (100s, 500s and 1000 s).

Step 8. Repeat steps 3-5 for (Room temperature, 50% - 90% intensity, 500 s).


Step 9. Repeat steps 3-5 for (Room temperature, 50% - 90% intensity, 100 s).

Step 10. Repeat steps 3-5 for (Room temperature, 0% - 50% intensity, 500 s).

Step 11. Repeat steps 3-5 for (Room temperature, 0% - 50% intensity, 100 s).

#### Analysis:

1. Check the shape of the LM-OSL curves for various parameter equations
2. For each case estimate the  $I_{max}$  and  $t_{max}$  and plot them versus the parameter that changes value.
3. For each case calculate the photo-ionization cross section of the main OSL component.
4. Check the impact of the LM-OSL measurements to the TL glow curve.



**Fig. 1:** Experimental parameters and dependent variables in LM-OSL measurements.

$$\sigma = \frac{2P}{c_f \cdot I_0 \cdot u_m^2 (\beta + 1)}. \quad (8)$$

For the stimulation light intensity expressed in W/cm<sup>2</sup> one obtains

$$c_f = \frac{\text{Joule}}{\text{s cm}^2} = \frac{6.24146 \times 10^{18} \text{ eV}}{\text{s cm}^2}. \quad (9)$$

For a given stimulation photon energy  $E_\phi$

$$c_f = \frac{6.24146 \times 10^{18}}{E_\phi} \left[ \frac{\text{eV}}{\text{eV s cm}^2} \right]. \quad (10)$$

Eventually

$$\sigma = \frac{2P \cdot E_\phi}{6.24146 \cdot 10^{18} \cdot u_m^2 I_0 \cdot (\beta + 1)}. \quad (11)$$

Total measurement time

Time corresponding to maximum intensity

Maximum stimulation intensity in W/cm<sup>2</sup>

**Fig. 2:** Methodology for the calculation of the photo-ionization cross section.