


PRACTICES

Quantity Determination of Antipyrine (Phenazone)



M.W.=188



1-phenyl-2,3-dimethyl-5-pyrazolone

100 mg of substance is dissolved in 18 ml of 10% sodium acetate. Therefore, the reaction turn back resulting from the strong reducer effect of HI will be prevent. 15 ml of 0.1 N I₂ solution is added and waited for 20 minutes in dark place by often shaking. 5 ml of CCl₄ solution is added. It is shaken to dissolve the precipitate. Excess of iod is titrated by 0.1 N Na₂S₂O₃ solution. The end point of the reaction is determined by the disappearance of the pink color in the CCl₄ layer.

1 ml 0.1 N I₂ is equivalent to 9.412 mg phenazone
 (15 x f I₂) - (b x f Na₂S₂O₃)..... A mg

A x 100/tartim= % phenazone

Reaction equation:



Standart solutions used in iodometric titration

0.1N I₂, 0.1N Na₂S₂O₃

Preperation of 0.1N I₂ Solution

0.1 N iod solution contains 12.69 g of iod in 1000 ml. 13 g of iod and 30 g of KI is dissolved in 40 ml of water and complete volume to 1000 ml with water.

Adjustment of 0.1N I₂ Solution to Arsentrioxide

In the erlenmeyer, 0.2 g of dried As₂O₃ is dissolved in 2 ml of 20% NaOH by heating if necessary. 40 ml of water and 0.2 ml of methyl orange is added to the solution. Until the colour turn red from yellow, dilue HCl is added. Then, the solution is neutralized by NaHCO₃. In order to prevent the

formation of HI acid, 2 g of weak base Na_2CO_3 is added. The primary standard is titrated by 0.1 N I_2 solution in the presence of 10 drops of starch solution until forming blue colour.

1 ml of 0.1N I_2is equivalent to 0.004946 g As_2O_3

Reaction equation:



Preparation and adjustment of 0.1N $\text{Na}_2\text{S}_2\text{O}_3$ Solution

$\text{Na}_2\text{S}_2\text{O}_3 \cdot 5\text{H}_2\text{O}$: 248.19

1000 ml of 0.1N $\text{Na}_2\text{S}_2\text{O}_3$ solution contains 24.82 g of $\text{Na}_2\text{S}_2\text{O}_3$.

26 g of $\text{Na}_2\text{S}_2\text{O}_3$ and 200 mg of Na_2CO_3 are dissolved in newly boiled and chilled distilled water and complete volume to 1000 ml.

In the adjustment, the reaction between thiosulfate ion and elemental iodine is utilized. Therefore, potassium iodate, potassium bromate, potassium bichromate, copper which release iodine or iodine from KI in acidic solution are used as standard.

Since iodine and thiosulfate solutions are used together in practice, firstly the factor of the iodine solution is determined. Then, the iodine solution is used as standard solution for determining the factor of thiosulfate solution.

10 ml of 0.1 N I_2 solution is titrated by 0.1 N $\text{Na}_2\text{S}_2\text{O}_3$ solution in the presence of 2 ml of starch solution until the blue colour disappears.

$V_1 \times f_1 = V_2 \times f_2$ equation is used.