

Experiment – 5 Airflow System –Jet Investigation Experiment

Aim of this Experiment

Jet Investigation duct has been designed for operation with the Airflow System. The duct allows students to investigate the phenomenon of jet attachment to an adjacent wall and the so called Coanda effect.

Experimental Set – up

Experimental system has some essential items needed for the experimental use. It features a large capacity **airflow system**, a **plenum chamber**, **multi-tube monometer** and **Bernoulli investigation duct**.

The **Airflow System** has been specifically designed to allow students to investigate a wide range of and low speed air flow phenomena and fundamental aerodynamics. Airflow System base unit consists of a large capacity variable speed centrifugal fan with a separate aerodynamically designed plenum chamber containing multiple screens, flow straightener and acceleration section. The fan and plenum chamber are connected by a length of flexible hose and this allows the two components to be arranged in a variety of convenient locations either at bench or floor level.

A large number of optional ducts may be attached to the **plenum discharge** that allow investigation of airflow on the positive side of the fan.

In addition there are additional optional items that attach to the suction or intake side of the fan.

The ability to utilise both the intake and discharge sides of the fan, together with a continuously expanding range of optional accessories makes the **Airflow System** a very flexible and cost effective unit.




General experimental setup of Airflow System: 1. Fan; 2. Fan Speed Control; 3. RCCB & MCB Box; 4. Fan Outlet; 5. Fan Inlet; 6. RCCB; 7. MCB




General experimental setup of Plenum Chamber: 8. Plenum Chamber; 9. Plenum Discharge; 10. Plenum Inlet

Multi-tube Manometer has been designed for operation with the Airflow System. However as a 16 tube manometer it may equally be used in any application that is within its pressure range. Last apparatus for the experimental setup is Bernoulli investigation duct.




General experimental setup of Multi-Tube Manometer: 1. Manometer Tube Couplings; 2. Side Clamps; 3. Angle Indicator; 4. Rubber Feet; 5. Manometer Tubes; 6. Reservoir Tapping; 7. Reservoir; 8. Reservoir Clamp; 9. Reservoir Track; 10. Marker



General experimental setup of Jet Attachment: 1. Splitter Block; 2. Locking Nut; 3. Attachment wall; 4. Jet nozzle; 5. Flexible Seal; 6. Wall extension; 7. Splitter Lock Screw; 8. Jet nozzle lock screw; 9. Mounting nuts; 10. Jet effective width “w”; 11. Wall extension lock screw

Schematic Diagram of inner radius, outer radius and radial pressure tapping centre lines



Theory

In a variety of aerodynamic experiments, such as flow around a cylinder and the flow around a wing the phenomena of separation of flow is commonly observed. The jet attachment experiment deals with a phenomenon, that is effectively, the reverse of separation, in that a jet, emerging close to a wall will attach to the wall and follow the contour.




Figure Jet Flow: 1. strong recirculating turbulent field; 2. effect of opening an aperture in the separation field

It is found that even if the wall curves appreciably the jet will still attach to its surface and follow the curve. Use of this phenomena was proposed in the 1930s by Henry Coanda, who made several devices that use this form Jet deflection. For this reason, the phenomena are often referred to as the Coanda effect.

It is known that a jet emerging into a field will entrain fluid adjacent to the point that it exits the orifice. In most cases, the jetty emerging is very turbulent and results in very strong mixing. In figure 6.1, the jet leaves the slot, and in an area adjacent the jet separates, and there is a strong recirculating turbulent field. The presence of the wall restricts mixing on this side and entrainment. This, together with the recirculation in the separated area and fluid being entrained to the right of the jet, both pull and push the jet back towards the wall.

The pressure in the separation field is lower than the atmospheric pressure on the opposite side of the jet. Figure 6.2 shows the effect of opening an aperture in the separation field. Due to the low pressure area is entrained from outside, and this has the effect of extending the reattachment point further down the wall. If the entrainment rate is increased, the reattachment point will move further and further down the wall.

Given sufficient entrainment the jet can effectively be detached permanently from the wall. Once attached to the wall as in Figure the jet is unlikely to detach as to do so would require the addition of fluid between the jet and the wall. Separation is unlikely then unless the curvature of the wall is large or there is an adverse (increasing) pressure gradient.

The description overleaf assumes two dimensional flow. If, however the jet is round for example then air entrained around the ends of the jet will also be entrained and cause separation sooner.

Due to the importance of these effects, the ratio of jet width to jet thickness is important when investigating this phenomenon.




Figure Aspect Ratio; $AR = (d/w)$

The “Coanda effect” can be used in fluid logic controls to act as a “flip-flop” control. These are analogous to electronic flip-flop switches that effectively stay in a set state unless an external signal acts upon them, at which point they switch over to the opposite state.




Figure The Mechanism of Coanda Effect

The jet emerging at the top of the diagram causes turbulent recirculation and entrainment on both sides as it emerges. If both control flap A and B are closed, then there is no preferred route for the flow. If control flap B is opened, then air will be entrained on this side and the jet will move towards the wall on Outlet A and will attach to the wall. Even if control flap B is closed the jet will remain attached to the wall if the jet and divider between outlet A and B have been designed correctly.

Only if control flap A is opened will the flow be diverted to outlet B. Again unless acted upon by control B, the flow will remain attached to wall B and exit through outlet B.

By using this phenomena, fluid switches like these can be constructed and developed into and/or gates to build up logic circuits.

Jet Attachment Investigation experiment allows students to investigate the phenomenon in detail and operate a flip-flop in this manner.

Experiment -5.1

Investigation of Jet Attachment

Procedure

It is necessary to connect the manometer to the **Plenum chamber** for reference purposes. The **airflow fans** should be started and the **Multi-Tube Manometer** levels should be adjusted.

Connection to the Airflow System


Care must be exercised when connecting the manometer to the airflow system. The following method is suggested to prevent the manometer liquid from being either driven out of the manometer tubes or drawn into the tubes connected to the accessories.

Before starting the fan, connect the pressure hoses to the accessory in use and to the manometer. Note that the two outer tubes (left and right) are not normally connected/used.

Set the manometer to the vertical or inclined condition as required and adjust the reservoir to about mid-height. Record the atmospheric datum or zero level. Then start the fan and slowly increase the speed, at the same time monitoring the manometer levels. As the pressures in the various tubes move up and/or down adjust the **reservoir** level also up or down, so that the liquid levels are kept within the range of the manometer.

Once the fan is running at the desired speed make any final adjustments to the reservoir level to set the atmospheric datum to a convenient value using the two outer tubes as a reference. Record this atmospheric datum as the reference value. It is this value that will be either taken from, or added to the other levels recorded on the manometer tubes.


Adjustment Details for the Jet Attachment Investigation accessory: Please also refer to the figure.




The width of the jet w and the offset of the jet y from the moveable wall can all be adjusted by moving the two steel **jet nozzle** plates. The plates are secured in position by the **jet nozzle lock screws**.

In order to examine jet attachment without the complication of the splitter and second wall the **splitterblock** can be removed by taking out the **locking screws** and sliding the block out of the duct.

Next one of the **attachment walls** may be swung to the horizontal position (as in the picture below on the right) to remove it from the flow.



In the above arrangement the jet has been located so that y , the distance between the edge of the jet and the left hand wall is 0mm. This can be checked with the angle gauge supplied as shown below.




Using a combination of the angle gauge and a rule, the distance y may be measured and adjusted. Note that the value of y may be set at positive or negative values. In the above diagram, a negative setting would be regarded as the jet edge being closer to the centre line of the wall pivot (as indicated by the arrows).

The **flexible seal** may either be allowed to sit against the outside of the wall or held back so that the gap between the jet and the wall is open. Tests are carried out with and without the seal in position in order to see how the induced flow referred to Figure 6 affects the turning angle of the air.

The procedure is to set the wall angle β to zero and bring the fan to a convenient speed close to maximum. This is then held constant through the tests by maintaining the plenum chamber pressure at a constant value.


Using a flat surface aligned with the flow, now slowly move the wall to increase the wall angle β .



As shown a rule is a convenient tool. By holding a hand adjacent to the end of the wall (at A above) the jet will be felt as it turns and follows the wall. As the wall angle β is increased a point will be found where the flow detaches and continues along arrow B. It may be necessary to repeat this process several times in order to determine the angle at which the flow “just breaks away” this angle is β_s .

The point at which the flow detaches can often also be “heard” as a change in the sound of the flow. By using the angle gauge supplied, measure the wall angle β_s at this condition.

Then slowly move the wall back to reduce the angle and determine the point at which the flow reattaches to the wall. This angle is β_r . This procedure should be repeated several times until the angle is determined and repeatable, to within about 1° .



It is recommended that to save setting time, the procedure is repeated with the **flexible seal** in position and not in position. It will be found that the results are very different in the two conditions.

The procedure is repeated at different valves of y offset and constant jet width w .

If time permits the whole procedure may be repeated at different jet width w values.

Typical Results


Jet Width $w = 10\text{mm}$

Aspect Ratio


$$\begin{aligned}
 AR &= \frac{d}{b} \\
 &= \frac{100}{10} \\
 &= 10
 \end{aligned}$$

y mm	Flexible seal Closed		Flexible Seal Open	
	β_s	β_r	β_s	β_r
-4	48	38	45	35
-2	50	40	50	38
0	53	43	53	43
2	64	47	51	34
4	66	49	36	25
6	71	56	30	23
8	73	58	27	20
10	77	30	21	18
14	73	25	17	14
16	67	25	16	13

The data is presented graphically



Data from more general test




Appendix – I Some Useful Data

Manometer Fluid Density 800 kg/m^3



Jet Depth	d	100mm
Jet width	w	Adjustable
Distance of jet from wall	x	12.0mm
Wall radius		R 6.0mm
Wall angle		β Adjustable
Edge of Jet to Wall distance	y	Adjustable



Appendix –II Symbols and Units

<u>UNITS</u>	<u>Symbol</u>	<u>Designation</u>	<u>Unit</u>
	C	Constant	-
	C _p	Non dimensional Pressure coefficient	-
	g	Acceleration due to gravity	m/s ²
	h	Manometer liquid height	mm
	p	Static pressure	kN/m ²
	P	Total pressure	kN/m ²
	ρ	Density of manometer fluid	Kg/m ³
	\dot{Q}	Flow Volume	m ³ /s
	r	radius dimension	m
	v	Local velocity	m/s
	V	Free Stream Velocity	m/s